

Human PBLD/MAWBP (His tag) recombinant protein



Catalog Number: 518455

General Information

Gene Name Synonym

Phenazine Biosynthesis-Like Domain-Containing Protein, MAWD-Binding Protein, MAWDBP, Unknown Protein 32 From 2D-PAGE of Liver Tissue, PBLD, MAWBP

Protein Construction

A DNA sequence encoding the Human PBLD/MAWBP Met1-Ala288 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 100mM NaCl, 1mM DTT, 30% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human PBLD/MAWBP (His tag) recombinant protein