

Human CCL3L1/D17S1718/G0S19-2 (His tag) recombinant protein



Catalog Number: 518461

General Information

Gene Name Synonym

C-C Motif Chemokine 3-Like 1, G0/G1 Switch
Regulatory Protein 19-2, LD78-Beta(1-70), PAT
464.2, Small-Inducible Cytokine A3-Like 1,
Tonsillar Lymphocyte LD78 Beta Protein, CCL3L1,
D17S1718, G0S19-2, SCYA3L1, CCL3L3

Protein Construction

A DNA sequence encoding the Human
CCL3L1/D17S1718/G0S19-2 Ala24-Ala93 is
expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL
method.

Stability

Store for up to 12 months at -20°C to -80°C as
lyophilized powder. Aliquots of reconstituted
samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 µm filtered solution of
20mM PB, 150mM NaCl, pH 7.2.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C,
though stable at room temperature for 3 weeks.
Reconstituted protein solution can be stored at
4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not
mix by vortex or pipetting. It is not recommended
to reconstitute to a concentration less than
100µg/ml. Dissolve the lyophilized protein in
ddH₂O. Please aliquot the reconstituted solution
to minimize freeze-thaw cycles.

SDS-PAGE



Human CCL3L1/D17S1718/G0S19-2 (His tag)
recombinant protein