

Human SAA4/CSAA (His tag) recombinant protein



Catalog Number: 518464

General Information

Gene Name Synonym

Serum Amyloid A-4 Protein, Constitutively Expressed Serum Amyloid A Protein, C-SAA, SAA4, CSAA

Protein Construction

A DNA sequence encoding the Human SAA4/CSAA Glu19-Thr130 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris, 200mM NaCl, 0.5mM EDTA, 20% Glycerol, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human SAA4/CSAA (His tag) recombinant protein