

Human ZNRF1/NIN283 (His tag) recombinant protein



Catalog Number: 518476

General Information

Gene Name Synonym

E3 ubiquitin-protein ligase ZNRF1, nerve injury-induced gene 283 protein, zinc/ring finger protein 1, NIN283, ZNRF1

Protein Construction

A DNA sequence encoding the Human ZNRF1/NIN283 Pro74-Thr178 is expressed with a 6His tag at the N-terminus, 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 0.15M NaCl, 1mMEDTA, 10% glycerol, pH 8.0 .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ZNRF1/NIN283 (His tag) recombinant protein