

Mouse Vegfd/Figf (His tag) recombinant protein



Catalog Number: 518490

General Information

Gene Name Synonym:

Vascular endothelial growth factor D, c-Fos-induced growth factor, FIGF, VEGFD,

Protein Construction:

A DNA sequence encoding the Mouse Vegfd/Figf Phe98-Ser206 is expressed with a 6His tag at the C-terminus.

Source: Mouse

Expression Host: Human Cells

QC Testing

Activity:

Purity:

>95% as determined by SDS-PAGE.

Endotoxin:

<1.0 EU per μg protein as determined by the LAL method.

Stability:

Store for up to 6 months at -20°C

Predicted N terminal:

Molecular Mass:

Formulation:

Supplied as a 0.2 μm filtered solution of 20mM PB, 150mM NaCl, pH7.4.


Usage Guide

Storage:

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

Reconstitution:

SDS-PAGE:



Mouse Vegfd/Figf (His tag) recombinant protein