

# Human IL22RA2/UNQ5793/PRO19598/PRO19822 (Fc tag) recombinant protein



Catalog Number: 518519

## General Information

### Gene Name Synonym

Interleukin-22 Receptor Subunit Alpha-2, IL-22 Receptor Subunit Alpha-2, IL-22R-Alpha-2, IL-22RA2, Cytokine Receptor Class-II Member 10, Cytokine Receptor Family 2 Member 10, CRF2-10, Cytokine Receptor Family Type 2 Soluble 1, CRF2-S1, Interleukin-22-Binding Protein, IL-22BP, IL22BP, ZcytoR16, IL22RA2

### Protein Construction

A DNA sequence encoding the Human IL22RA2/UNQ5793/PRO19598/PRO19822 Thr22-Pro231 is expressed with a Fc tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

### Stability

Store for up to 12 months at -20°C to -80°C as

lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

### Formulation

Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100µg/ml. Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

### SDS-PAGE



Human IL22RA2/UNQ5793/PRO19598/PRO19822 (Fc tag) recombinant protein