

Human HPGD/PGDH1/SDR36C1 (His tag) recombinant protein



Catalog Number: 518533

General Information

Gene Name Synonym

15-Hydroxyprostaglandin Dehydrogenase
[NAD(+)], 15-PGDH, Prostaglandin
Dehydrogenase 1, HPGD, PGDH1

Protein Construction

A DNA sequence encoding the Human
HPGD/PGDH1/SDR36C1 Met1-Gln266 is
expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM
HEPES, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human HPGD/PGDH1/SDR36C1 (His tag)
recombinant protein