

# Human TRAT1/TCRIM/HSPC062 (His tag) recombinant protein



Catalog Number: 518647

## General Information

### Gene Name Synonym

T-Cell Receptor-Associated Transmembrane Adapter 1, T-Cell Receptor-Interacting Molecule, TRIM, pp29/30, TRAT1, TCRIM, HSPC062

### Protein Construction

A DNA sequence encoding the Human TRAT1/TCRIM/HSPC062 Asn29-Asn186 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human TRAT1/TCRIM/HSPC062 (His tag)  
recombinant protein