

Human FABP6/ILBP/ILLBP (His tag) recombinant protein



Catalog Number: 518655

General Information

Gene Name Synonym

Gastrotropin, GT, Fatty Acid-Binding Protein 6, Ileal Lipid-Binding Protein, ILBP, Intestinal 15 kDa Protein, I-15P, Intestinal Bile Acid-Binding Protein, I-BABP, FABP6, ILBP, ILLBP

Protein Construction

A DNA sequence encoding the Human FABP6/ILBP/ILLBP Met1-Ala128 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 0.5mM DTT, 50% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human FABP6/ILBP/ILLBP (His tag) recombinant protein