

# Human SH3KBP1/CIN85 (His tag) recombinant protein



Catalog Number: 518667

## General Information

### Gene Name Synonym

SH3 Domain-Containing Kinase-Binding Protein 1, CD2-Binding Protein 3, CD2BP3, Cbl-Interacting Protein of 85 kDa, Human Src Family Kinase-Binding Protein 1, HSB-1, SH3KBP1, CIN85

### Protein Construction

A DNA sequence encoding the Human SH3KBP1/CIN85 Met1-Lys665 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted samples are stable at  $-20^{\circ}\text{C}$  for 3 months.

## Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of 50mM Tris-HCl, 10mM reduced Glutathione, pH=8.0.

## Usage Guide


### Storage

Lyophilized protein should be stored at  $<-20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human SH3KBP1/CIN85 (His tag) recombinant protein