

# Human SENP8/DEN1/NEDP1/PRSC2/FKSG8 (His tag) recombinant protein



Catalog Number: 518676

## General Information

### Gene Name Synonym

Sentrin-Specific Protease 8, Deneddylase-1, NEDD8-Specific Protease 1, Protease Cysteine 2, Sentrin/SUMO-Specific Protease SENP8, SENP8, DEN1, NEDP1, PRSC2

### Protein Construction

A DNA sequence encoding the Human SENP8/DEN1/NEDP1/PRSC2/FKSG8 Met1-Lys212 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM PB, 150mM NaCl, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human SENP8/DEN1/NEDP1/PRSC2/FKSG8 (His tag) recombinant protein