

Human SMAD2/MADH2/MADR2 (His & Flag tag) recombinant protein



Catalog Number: 518683

General Information

Gene Name Synonym

Mothers against decapentaplegic homolog 2; MAD homolog 2; Mad-related protein 2; hMAD-2; SMAD family member 2; SMAD 2; Smad2; Hsmad2; SMAD2; MADH2; MADR2

Protein Construction

A DNA sequence encoding the Human SMAD2/MADH2/MADR2 Ser2-Ser467 is expressed with a 6His, Flag tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at $< -20^{\circ}\text{C}$ for 3 months.

Formulation

Lyophilized from a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 500mM NaCl, pH7.5.

Usage Guide


Storage

Lyophilized protein should be stored at $< -20^{\circ}\text{C}$, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at $4-7^{\circ}\text{C}$ for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than $100\ \mu\text{g}/\text{ml}$. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human SMAD2/MADH2/MADR2 (His & Flag tag)
recombinant protein