

Human HYAL1/LUCA1 (His tag) recombinant protein



Catalog Number: 518691

General Information

Gene Name Synonym

Hyaluronidase-1, Hyal-1,
Hyaluronoglucosaminidase-1, Lung Carcinoma
Protein 1, LuCa-1, HYAL1, LUCA1

Protein Construction

A DNA sequence encoding the Human
HYAL1/LUCA1 Phe22-Trp435 is expressed with a
6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM
TrisHCl, 150mM NaCl, 10%Glycerol, pH7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human HYAL1/LUCA1 (His tag) recombinant
protein