

Human PCSK9/NARC1/PSEC0052 (His tag) recombinant protein



Catalog Number: 518702

General Information

Gene Name Synonym

Proprotein Convertase Subtilisin/Kexin Type 9,
Neural Apoptosis-Regulated Convertase 1,
NARC-1, Proprotein Convertase 9, PC9,
Subtilisin/Kexin-Like Protease PC9, PCSK9,
NARC1

Protein Construction

A DNA sequence encoding the Human
PCSK9/NARC1/PSEC0052 Gln31-Gln692 is
expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 50mM
HEPES, 150mM NaCl, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PCSK9/NARC1/PSEC0052 (His tag)
recombinant protein