

# Human FLRT2/KIAA0405/UNQ232/PRO265 (His tag) recombinant protein



Catalog Number: 518704

## General Information

### Gene Name Synonym

Leucine-Rich Repeat Transmembrane Protein  
FLRT2, Fibronectin-Like Domain-Containing  
Leucine-Rich Transmembrane Protein 2, FLRT2,  
KIAA0405

### Protein Construction

A DNA sequence encoding the Human  
FLRT2/KIAA0405/UNQ232/PRO265 Cys36-Ser539  
is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL  
method.

### Stability

Store for up to 12 months at -20°C to -80°C as  
lyophilized powder. Aliquots of reconstituted  
samples are stable at < -20°C for 3 months.

## Formulation

Lyophilized from a 0.2  $\mu$ m filtered solution of  
20mM PB, 150mM NaCl, pH 7.2.

## Usage Guide


### Storage

Lyophilized protein should be stored at <-20°C,  
though stable at room temperature for 3 weeks.  
Reconstituted protein solution can be stored at  
4-7°C for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not  
mix by vortex or pipetting. It is not recommended  
to reconstitute to a concentration less than  
100 $\mu$ g/ml. Dissolve the lyophilized protein in  
ddH<sub>2</sub>O. Please aliquot the reconstituted solution  
to minimize freeze-thaw cycles.

## SDS-PAGE



Human FLRT2/KIAA0405/UNQ232/PRO265 (His  
tag) recombinant protein