

General Information

Gene Name Synonym

E3 Ubiquitin-Protein Ligase CHIP, Antigen NY-CO-7, CLL-Associated Antigen KW-8, Carboxy Terminus of Hsp70-Interacting Protein, STIP1 Homology and U Box-Containing Protein 1, STUB1, CHIP

Protein Construction

A DNA sequence encoding the Human STUB1/CHIP/PP1131 Met1-Tyr303 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human STUB1/CHIP/PP1131 recombinant protein