

Human BID recombinant protein



Catalog Number: 518719

General Information

Gene Name Synonym

BH3-Interacting Domain Death Agonist, p22 BID, BID

Protein Construction

A DNA sequence encoding the Human BID Met1-Asp195 is expressed.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 100mM KCl, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human BID recombinant protein