

Human PPIE/CYP33 (His tag) recombinant protein



Catalog Number: 518738

General Information

Gene Name Synonym

Peptidyl-Prolyl Cis-Trans Isomerase E, PPIase E, Cyclophilin E, Cyclophilin-33, Rotamase E, PPIE, CYP33

Protein Construction

A DNA sequence encoding the Human PPIE/CYP33 Met1-Val301 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human PPIE/CYP33 (His tag) recombinant protein