

Human GAS7/KIAA0394 (His tag) recombinant protein



Catalog Number: 518742

General Information

Gene Name Synonym

Growth Arrest-Specific Protein 7, GAS-7, GAS7, KIAA0394

Protein Construction

A DNA sequence encoding the Human GAS7/KIAA0394 Met1-Ile412 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 100mM NaCl, 2mM DTT, 10% Glycerol, pH 8.8.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human GAS7/KIAA0394 (His tag) recombinant protein