

Catalog Number: 518769

## General Information

### Gene Name Synonym

Ubiquitin-like-conjugating enzyme ATG3;  
Autophagy-related protein 3; APG3-like; hApg3;  
Protein PC3-96; ATG3; APG3; APG3L

### Protein Construction

A DNA sequence encoding the Human  
ATG3/APG3/APG3L Met1-Met314 is expressed.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM  
TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.0.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human ATG3/APG3/APG3L recombinant protein