

Human CDK2AP1/CDKAP1/DOC1 (His tag) recombinant protein



Catalog Number: 518781

General Information

Gene Name Synonym

Cyclin-Dependent Kinase 2-Associated Protein 1, CDK2-Associated Protein 1, Deleted in Oral Cancer 1, DOC-1, Putative Oral Cancer Suppressor, CDK2AP1, CDKAP1, DOC1

Protein Construction

A DNA sequence encoding the Human CDK2AP1/CDKAP1/DOC1 Met1-Ser115 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 12 months at -20°C to -80°C as lyophilized powder. Aliquots of reconstituted samples are stable at < -20°C for 3 months.

Formulation

Lyophilized from a 0.2 μ m filtered solution of 20mM Tris, 150mM NaCl, pH8.0.

Usage Guide


Storage

Lyophilized protein should be stored at <-20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days.

Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than 100 μ g/ml. Dissolve the lyophilized protein in ddH₂O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

SDS-PAGE



Human CDK2AP1/CDKAP1/DOC1 (His tag)
recombinant protein