

# Human GNMT (His tag) recombinant protein



Catalog Number: 518794

## General Information

### Gene Name Synonym

Glycine N-Methyltransferase, GNMT

### Protein Construction

A DNA sequence encoding the Human GNMT Met1-Asg294 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E.coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human GNMT (His tag) recombinant protein