

# Human HMOX1/HO/HO1 recombinant protein



Catalog Number: 518799

## General Information

### Gene Name Synonym

Heme Oxygenase 1, HO-1, HMOX1, HO, HO1

### Protein Construction

A DNA sequence encoding the Human HMOX1/HO/HO1 Met1-Thr261 is expressed.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

## Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

## Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human HMOX1/HO/HO1 recombinant protein