

# Human TXNDC12/TLP19/UNQ713/PRO1376 (His tag) recombinant protein



Catalog Number: 518885

## General Information

### Gene Name Synonym

Thioredoxin Domain-Containing Protein 12, Endoplasmic Reticulum Resident Protein 18, ER Protein 18, ERp18, Endoplasmic Reticulum Resident Protein 19, ER Protein 19, ERp19, Thioredoxin-Like Protein p19, hTLP19, TXNDC12, TLP19

### Protein Construction

A DNA sequence encoding the Human TXNDC12/TLP19/UNQ713/PRO1376 His27-Leu168 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human TXNDC12/TLP19/UNQ713/PRO1376 (His tag) recombinant protein