

# Human VCP (His tag) recombinant protein



Catalog Number: 518930

## General Information

### Gene Name Synonym

Transitional Endoplasmic Reticulum ATPase, TER ATPase, 15S Mg(2+)-ATPase p97 Subunit, Valosin-Containing Protein, VCP

### Protein Construction

A DNA sequence encoding the Human VCP Ala2-Asn589 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 50mM TrisHCl pH8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human VCP (His tag) recombinant protein