

# Cavia porcellus CTLA4 (His tag) recombinant protein



Catalog Number: 518934

## General Information

### Gene Name Synonym

Cytotoxic T-lymphocyte protein 4; Cytotoxic T-lymphocyte-associated antigen 4; CTLA-4; CD152; CTLA4

### Protein Construction

A DNA sequence encoding the Cavia porcellus CTLA4 Ala37-Asp161 is expressed with a 6His tag at the C-terminus.

### Organism

Cavia porcellus

### Expression Host

Yeast

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of PBS, pH7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Cavia porcellus CTLA4 (His tag) recombinant protein