

# Mouse IL-27 (His tag) recombinant protein



Catalog Number: 518984

## General Information

### Gene Name Synonym

Interleukin-27 subunit alpha,IL-27 subunit alpha,IL-27-A,IL27-A,p28,Il27,Il27a,Interleukin-27 subunit beta,Ebi3,IL-27 subunit beta,IL-27B,Epstein-Barr virus-induced gene 3 protein homolog,Ebi3,Il27b

### Protein Construction

A DNA sequence encoding the Mouse IL-27 Tyr19-Pro228&Phe29-Ser234 is expressed

### Organism

Mouse

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

## Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of PBS,pH7.4.

## Usage Guide


### Storage

Lyophilized protein should be stored at  $<-20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Mouse IL-27 (His tag) recombinant protein