

Human SULT2B1/HSST2 (His tag) recombinant protein



Catalog Number: 519006

General Information

Gene Name Synonym

Sulfotransferase Family Cytosolic 2B Member 1, ST2B1, Sulfotransferase 2B1, Alcohol Sulfotransferase, Hydroxysteroid Sulfotransferase 2, SULT2B1, HSST2

Protein Construction

A DNA sequence encoding the Human SULT2B1/HSST2 Met1-Glu311 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 0.1M NaCl, 1mM EDTA, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human SULT2B1/HSST2 (His tag) recombinant protein