

Human UBE2A/RAD6A (GST & His tag) recombinant protein



Catalog Number: 519021

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 A, RAD6 Homolog A, HR6A, hHR6A, Ubiquitin Carrier Protein A, Ubiquitin-Protein Ligase A, UBE2A, RAD6A

Protein Construction

A DNA sequence encoding the Human UBE2A/RAD6A Met1-Cys152 is expressed with a GST, 6His tag at the N-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2A/RAD6A (GST & His tag)
recombinant protein