

Human LDHB (His tag) recombinant protein



Catalog Number: 519029

General Information

Gene Name Synonym

L-lactate Dehydrogenase B Chain, LDH-B, LDH Heart Subunit, LDH-H, Renal Carcinoma Antigen NY-REN-46, LDHB

Protein Construction

A DNA sequence encoding the Human LDHB Met1-Leu334 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 100mM NaCl, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human LDHB (His tag) recombinant protein