

# Human TPP1/CLN2/GIG1/UNQ267/PRO304 (His tag) recombinant protein



Catalog Number: 519087

## General Information

### Gene Name Synonym

Tripeptidyl-Peptidase 1, TPP-1, Cell Growth-Inhibiting Gene 1 Protein, Lysosomal Pepstatin-Insensitive Protease, LPIC, Tripeptidyl Aminopeptidase, Tripeptidyl-Peptidase I, TPP-I, TPP1, CLN2

### Protein Construction

A DNA sequence encoding the Human TPP1/CLN2/GIG1/UNQ267/PRO304 Ser20-Pro563 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM GaCl<sub>2</sub>, 10% Glycerol, pH7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human TPP1/CLN2/GIG1/UNQ267/PRO304 (His tag) recombinant protein