

Human SELENOI/EPT1/KIAA1724/SELI (GST tag) recombinant protein



Catalog Number: 519105

General Information

Gene Name Synonym

Ethanolaminephosphotransferase 1, hEPT1, Selenoprotein I, SelI, EPT1, KIAA1724, SELI

Protein Construction

A DNA sequence encoding the Human SELENOI/EPT1/KIAA1724/SELI Met1-Pro50 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human SELENOI/EPT1/KIAA1724/SELI (GST tag)
recombinant protein