

Mouse Timp2/Timp-2 (His tag) recombinant protein



Catalog Number: 519113

General Information

Gene Name Synonym

TIMP-2; CSC-21K tissue inhibitor of metalloproteinase 2; metalloproteinase inhibitor 2; TIMP metalloproteinase inhibitor 2; Tissue inhibitor of metalloproteinase 2.

Protein Construction

A DNA sequence encoding the Mouse Timp2/Timp-2 Cys27-Pro220 is expressed

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Timp2/Timp-2 (His tag) recombinant protein