

General Information

Gene Name Synonym

Alpha-2-antiplasmin, Alpha-2-AP, Alpha-2-plasmin inhibitor, Alpha-2-PI, Serpin F2, Serpinf2

Protein Construction

A DNA sequence encoding the Mouse Serpinf2/Pli Val28-Lys491 is expressed with a Fc tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM Tris, 150mM NaCl, pH 8.0 (20% glycerol, 1mM DTT).

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Serpinf2/Pli (Fc tag) recombinant protein