

# Human UBE2T/HSPC150/PIG50 (His tag) recombinant protein



Catalog Number: 519242

## General Information

### Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 T, Cell Proliferation-Inducing Gene 50 Protein, Ubiquitin Carrier Protein T, Ubiquitin-Protein Ligase T, UBE2T

### Protein Construction

A DNA sequence encoding the Human UBE2T/HSPC150/PIG50 Met1-Val197 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human UBE2T/HSPC150/PIG50 (His tag)  
recombinant protein