

Human ZNF100 (His tag) recombinant protein



Catalog Number: 519255

General Information

Gene Name Synonym

Zinc Finger Protein 100, ZNF100

Protein Construction

A DNA sequence encoding the Human ZNF100 Met1-Lys109 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 200mM NaCl, 50mM Imidazole, 1mM ZnCl₂, 30% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ZNF100 (His tag) recombinant protein