

# Human ODC1 (His tag) recombinant protein



Catalog Number: 519286

## General Information

### Gene Name Synonym

Ornithine decarboxylase, ODC1,

### Protein Construction

A DNA sequence encoding the Human ODC1 Met1-Val461 is expressed with a T7 tag at the N-terminus, 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

## Stability

Store for up to 6 months at -20°C

## Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 8.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human ODC1 (His tag) recombinant protein