

# Human POMGNT1/MGAT1.2/UNQ746/PRO1475 (His tag) recombinant protein



Catalog Number: 519310

## General Information

### Gene Name Synonym

Protein O-Linked-Mannose Beta-1,2-N-Acetylglucosaminyltransferase 1, POMGnT1, UDP-GlcNAc:Alpha-D-Mannoside Beta-1,2-N-Acetylglucosaminyltransferase I.2, GnT I.2, POMGNT1, MGAT1.2

### Protein Construction

A DNA sequence encoding the Human POMGNT1/MGAT1.2/UNQ746/PRO1475 Leu59-Thr660 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human POMGNT1/MGAT1.2/UNQ746/PRO1475 (His tag) recombinant protein