

Human B4GALT4/UNQ552/PRO1109 (His tag) recombinant protein



Catalog Number: 519311

General Information

Gene Name Synonym

Beta-1,4-Galactosyltransferase 4, Beta-1,4-GalTase 4, Beta4Gal-T4, b4Gal-T4, UDP-Gal:Beta-GlcNAc Beta-1,4-Galactosyltransferase 4, UDP-Galactose:Beta-N-Acetylglucosamine Beta-1,4-Galactosyltransferase 4, B4GALT4

Protein Construction

A DNA sequence encoding the Human B4GALT4/UNQ552/PRO1109 Gln39-Ala344 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 150mM NaCl, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human B4GALT4/UNQ552/PRO1109 (His tag)
recombinant protein