

# Human UBE2V2/MMS2/UEV2 (His tag) recombinant protein



Catalog Number: 519321

## General Information

### Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 Variant 2, DDVit 1, Enterocyte Differentiation-Associated Factor 1, EDAF-1, Enterocyte Differentiation-Promoting Factor 1, EDPF-1, MMS2 Homolog, Vitamin D3-Inducible Protein, UBE2V2, MMS2, UEV2

### Protein Construction

A DNA sequence encoding the Human UBE2V2/MMS2/UEV2 Met1-Asn145 is expressed with a 6His tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu$ g protein as determined by the LAL method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2  $\mu$ m filtered solution of 50mm HEPES, 150mM NaCl, pH 7.0.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human UBE2V2/MMS2/UEV2 (His tag)  
recombinant protein