

Mouse Thbd (His tag) recombinant protein



Catalog Number: 519328

General Information

Gene Name Synonym

Thrombomodulin; TM; Fetomodulin; CD141; BDCA-3; Thbd

Protein Construction

A DNA sequence encoding the Mouse Thbd Leu17-Ser517 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM Tris, 150mM NaCl, pH8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Thbd (His tag) recombinant protein