

Human TPH2 (His tag) recombinant protein



Catalog Number: 519350

General Information

Gene Name Synonym

Tryptophan 5-Hydroxylase 2, Neuronal
Tryptophan Hydroxylase, Tryptophan 5-
Monooxygenase 2, TPH2, NTPH

Protein Construction

A DNA sequence encoding the Human TPH2
Met1-Ile496 is expressed with a 6His tag at the N-
terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB \square
150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human TPH2 (His tag) recombinant protein