

Human TMPRSS11B/HATL5 (His tag) recombinant protein



Catalog Number: 519356

General Information

Gene Name Synonym

Transmembrane Protease Serine 11B, Airway
Trypsin-Like Protease 5, TMPRSS11B, HATL5

Protein Construction

A DNA sequence encoding the Human
TMPRSS11B/HATL5 Leu39-Leu416 is expressed
with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL
method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM
PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human TMPRSS11B/HATL5 (His tag) recombinant
protein