

Human NEU1/NANH (His tag) recombinant protein



Catalog Number: 519365

General Information

Gene Name Synonym

Sialidase-1, Acetylneuraminyl Hydrolase, G9 Sialidase, Lysosomal Sialidase, N-Acetyl-Alpha-Neuraminidase 1, NEU1, NANH

Protein Construction

A DNA sequence encoding the Human NEU1/NANH Glu48-Leu415 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human NEU1/NANH (His tag) recombinant protein