

# Human IFNAR2/IFNABR/IFNARB (His tag) recombinant protein



Catalog Number: 519375

## General Information

### Gene Name Synonym

Interferon Alpha/Beta Receptor 2, IFN-R-2, IFN-Alpha Binding Protein, IFN-Alpha/Beta Receptor 2, Interferon Alpha Binding Protein, Type I Interferon Receptor 2, IFNAR2, IFNABR, IFNARB

### Protein Construction

A DNA sequence encoding the Human IFNAR2/IFNABR/IFNARB Ile27-Lys243 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 12 months at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$  as lyophilized powder. Aliquots of reconstituted samples are stable at  $< -20^{\circ}\text{C}$  for 3 months.

## Formulation

Lyophilized from a  $0.2\ \mu\text{m}$  filtered solution of 20mM PB, 150mM NaCl, pH 7.2.

## Usage Guide


### Storage

Lyophilized protein should be stored at  $< -20^{\circ}\text{C}$ , though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at  $4-7^{\circ}\text{C}$  for 2-7 days.

### Reconstitution

Always centrifuge tubes before opening. Do not mix by vortex or pipetting. It is not recommended to reconstitute to a concentration less than  $100\ \mu\text{g}/\text{ml}$ . Dissolve the lyophilized protein in ddH<sub>2</sub>O. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.

## SDS-PAGE



Human IFNAR2/IFNABR/IFNARB (His tag)  
recombinant protein