

# Human STAT5B (His tag) recombinant protein



Catalog Number: 519413

## General Information

### Gene Name Synonym

Signal Transducer and Activator of Transcription 5B, STAT5B

### Protein Construction

A DNA sequence encoding the Human STAT5B Met1-Thr321 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL

method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of PBS, 50% Glycerol, 1mM DTT, pH 7.4.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human STAT5B (His tag) recombinant protein