

Human GALK2/GK2 (His tag) recombinant protein



Catalog Number: 519427

General Information

Gene Name Synonym

N-acetylgalactosamine kinase, also known as GalNAc kinase, Galactokinase 2, GALK2 and GK2, is a member of the GHMP kinase family and GalK subfamily.

Protein Construction

A DNA sequence encoding the Human GALK2/GK2 Met1-Ala458 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μ g protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human GALK2/GK2 (His tag) recombinant protein