

Human SEMA3A/SEMAD (His tag) recombinant protein



Catalog Number: 519440

General Information

Gene Name Synonym

Synonyms: Semaphorin III, SEMA3A, SEMAD, Hsema-I, SEMA1, SemD, Semad, coll-1 , Hsema-III, Sema III

Protein Construction

A DNA sequence encoding the Human SEMA3A/SEMAD Asn21-Val771 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris 150mM NaCl 1mM EDTA 20%Glycerol, pH8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human SEMA3A/SEMAD (His tag) recombinant protein