

Human BDH2/DHRS6/SDR15C1/UNQ6308 (His tag) recombinant protein



Catalog Number: 519468

General Information

Gene Name Synonym

3-Hydroxybutyrate Dehydrogenase Type 2, Dehydrogenase/Reductase SDR Family Member 6, Oxidoreductase UCPA, R-Beta-Hydroxybutyrate Dehydrogenase, BDH2, DHRS6

Protein Construction

A DNA sequence encoding the Human BDH2/DHRS6/SDR15C1/UNQ6308 Met1-Leu245 is expressed with a 6His tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 100mM NaCl, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human BDH2/DHRS6/SDR15C1/UNQ6308 (His tag) recombinant protein