

Human GLUL/GLNS (His tag) recombinant protein



Catalog Number: 519492

General Information

Gene Name Synonym

Glutamine Synthetase, GS, Glutamate Decarboxylase, Glutamate--Ammonia Ligase, GLUL, GLNS

Protein Construction

A DNA sequence encoding the Human GLUL/GLNS Thr2-Asn373 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E.coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 200mM NaCl, 50mM Imidazole, pH 8.0.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human GLUL/GLNS (His tag) recombinant protein